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REF. 20304

Reinforced neck brace

DESIGNED AND DEVELOPED IN SWITZERLAND

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EN – The descriptions and images in this document are for example and commercial purposes only. Orthoservice reserves the right to modify them according to its needs.

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90DF6001491HT03.0622

 Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland

Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89

info@orthoservice.com · www.orthoservice.com

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Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13

info@orthoservice.de · www.orthoservice.de

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Tel.: +48 (0) 34 340 13 10 · (NIP): 9492246785

www.orthoservice.pl · info@orthoservice.pl

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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Reinforced neck brace
DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medcal device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

Code	REF.20304			
Size	S	M	L	XL
Neck circumf. cm	30/35	36/41	42/47	48/53
Height cm	6-8-10-12	6-8-10-12	6-8-10-12	6-8-10-12
Colour	grey			

MAINTENANCE
☒ Do not bleach
☒ No chemical cleaning
☒ Do not iron
☒ Do not tumble-dry
☒ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Mild or moderate distortions of the cervical spine (whiplash).
- Neck-induced vertigo
- Cervical osteoarthritis
- Moderate discopathies with associated cervicobrachial pain syndrome.

CONTRAINDICATIONS
Currently no known

CHARACTERISTICS AND MATERIALS

- Core made of semi-rigid foam ((with rounded edges) reinforced with polyethylene inside the brace
- Tubular jersey cover
- Velcro® fastening at back of neck
- Fully washable
- X-ray transparent
- Anatomically shaped design.

PUTTING ON THE APPLIANCE
A Hold the neck brace in a position so that the concave part is at the top and the strap is on the outside of the brace once it has been fastened.
B Rest the concave part at the front of the neck brace under the chin.
C Fasten the neck brace at the back using the Velcro® strap provided.

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