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SINGLE PATIENT
MULTIPLE USE

PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Arm sling

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.8403
Size	universal
Colour	black

Fits right and left

MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry
- Washing instructions: hand wash at max. 30°C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

MATERIALS

Polyamide, polypropylene and polyester.

INDICATIONS

- For use as a support for plaster bandages
- After injuries
- After operations

CONTRAINDICATIONS

- Currently no known

CHARACTERISTICS AND MATERIALS

- Nylon strap, length adjustable using the appropriate buckles
- Neck padding made from double-sided foam with a soft jersey fabric cover

PUTTING ON THE APPLIANCE

- 1 Place the strap around the back of the neck and let the two ends hang down at the front (fig. A).
- 2 Slide the forearm through the first looped end as far as the elbow.
- 3 Slide the forearm into the second looped end so the wrist rests in it.
- 4 Adjust the length using the buckles provided (fig. B). The arm should be bent almost at a right angle.
- 5 Position the protective padding around the neck (fig. C and D).

IF IT IS PREFERABLE TO HAVE LESS PRESSURE ON THE NECK, THE SLING CAN ALSO BE WORN BY PASSING THE STRAP UNDER THE ARMPIT.

- 1 Place the strap around the back of the neck and let the end of the strap on the opposite side to the supported arm hang down at the front.
- 2 Pass the other end underneath the shoulder of the supported arm and slide in the arm as far as the elbow (fig. E)
- 3 Slide the arm into the remaining free end of the strap so the wrist rests in the loop (fig. F)
- 4 Adjust the length using the buckles provided (fig. G). The arm should be bent almost at a right angle.
- 5 Position the protective padding around the neck (fig. H)

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