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SINGLE PATIENT
MULTIPLE USE

PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Rigid lumbar orthosis

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

| | | | | |
|-------------------------|---------------|--------|---------|---------|
| Code | REF.8909 | | | |
| Size | S | M | L | XL |
| Pelvis circumference cm | 75/85 | 85/100 | 100/120 | 120/135 |
| Colour | beige / cream | | | |

MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry
- Washing instructions: wipe the plastic parts with a soft sponge dipped in lukewarm water and mild soap. Dry with a cloth. Hand wash the padding and the straps in tepid 30°C water with mild soap. Leave to dry away from direct heat.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Outcomes of lumbosacral vertebral fractures
- Spondylolysis and L4-L5 listesis
- Useful in the post-operative treatment of lumbar arthrosis

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Structure in KYDEX® plastic material
- Two sections: lumbosacral and abdominal
- Joining belts with Velcro® closure
- Soft coupled panels combined with the rigid structure

PUTTING ON THE APPLIANCE

- 1 If necessary, mould the "Kydex" plastic lumbar frame with the help of a heat gun.
- 2 Insert the straps into the abdominal panel.
- 3 Tighten the straps beginning with the central one, then the lower ones and then the upper ones.

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Velcro® est une marque enregistrée de Velcro Industries B.V.