ice

REF. 90000

Patella-Sehnenband
Patella strap
Bande rotulienne
Коленный ремень
Cinturino sottrotuleo

ORTHOSERVICE RO+TEN

Headquarter: ORTHOSERVICE AG
Via Milano 7 - CH-6830 Chiasso (TI) - Switzerland
Tel. 0041 91 8220088 - Fax 0041 91 8220089
info@orthoservice.com - www.orthoservice.com

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ENGLISH

Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

SELECTION/SIZE

Size: Universal

Colour: black

Fits right and left

MAINTENANCE

Do not bleach No chemical cleaning

Do not iron Do not tumble-dry

Washing instructions: hand wash at max. 30° C with neutral soap, do not dry in the vicinity of heat sources

Do not dispose of the product or any of its components into the environment.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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ENGLISH

FEATURES AND MATERIALS

- Elastic front strap for pressure adjustment. Ensures constant pressure, even when moving
- Rigid Velcro® surface for stable fixing of the pressure pad or the cooling pad (optional)
- Fused foam pressure pad with skin-friendly fabric cover. Imprinted markings for precise positioning. Imprinted markings for precise positioning. The skin-friendly gel element concentrates pressure on the tendon and protects against slippage
- Padded Velcro® strap with symmetrical tensioning system for easy positioning. Soft border for pleasant wearing comfort, especially at the back of the knee
- Optional: Pad for cooling therapy, with Velcro® attachment to exchange it easily with the pressure pad

INDICATIONS

- Tendinitis and tendinopathy in the patella region
- Morbus Osgood-Schlatter (inflammation of the front tibial apophysis)

CONTRAINDICATIONS

Currently no known

PUTTING ON THE APPLIANCE

- Pull one of the two ends of the strap out of the appropriate loop (fig. A).
- If the pelotte with the pressure pad has not already been positioned, affix this to the securing band using the adhesive points (fig. B).
- Put on the strap in such a way that the pressure pad is located below the knee at the painful position and the cavity of the gel pressure pad is facing upwards (fig. C).
- Feed the free end of the strap (the one you pulled out earlier) back through the appropriate loop (fig. D).
- Also unfasten the end of the Velcro® closure and pull the strap tight evenly (fig. E).
- Secure the ends to the strap using the Velcro® closure (fig. F).
- Check the position of the pressure pad.
- If the strap is uncomfortable, or causes a tingling sensation, loosen the closure a little. If the strap is too loose, it may be tightened up further.

COOLING

- Take the gel pad out of the velour pocket.
- Prior to use, put the gel pad in the freezer for at least two to three hours. (NOTE: DO NOT put the velour pocket in the freezer with the gel pad inserted).
- How to use:
 - Prior to use, put the gel pad into the velour pocket (fig. G). Place the gel pad and the pocket on the Velcro® surface in place of the pressure pad that is used for active usage (fig. H and I)
 - The strap should be applied in exactly the same way as described above
 - Application time: As prescribed by the physician

WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant that country's regulations concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

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