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1 Widen die Bandage (without pulling it out from the ring) and fit it (pull up from the hand towards the arm) (fig. C).

2 Place the bandage 3 or 4 fingers below the epicondyle (fig. D) and turn until the pressure pad is located on the extensor muscles of the hand (fig. E).

3 Open the Velcro® closure and tighten the bandage until the white line on the inside of the ring appears (connected to the elastic band) (fig. F).

4 Close the bandage with the Velcro® closure (fig. G).

5 Check that the pressure pad is in the correct position (point 3).

6 Should the bandage disturb too much when wearing it or cause a tingling sensation, slightly loosen the closure.

PUTTING ON THE APPLIANCE FOR THE PATIENT

1 Widen the bandage (without pulling it out from the ring) and fit it (pull up from the hand towards the arm) (fig. C).

2 Place the bandage 3 or 4 fingers below the epicondyle (fig. D) and turn until the pressure pad is located on the extensor muscles of the hand (fig. E).

3 Open the Velcro® closure and tighten the bandage until the white line on the inside of the ring appears (connected to the elastic band) (fig. F).

4 Close the bandage with the Velcro® closure (fig. G).

5 Check that the pressure pad is in the correct position (point 3).

6 Should the bandage disturb too much when wearing it or cause a tingling sensation, slightly loosen the closure.

COLDING

7 Take the gel cushion out of its velour cover.

8 Place the cushion into the freezer at least 2-3 hours prior to its application (NOTE: Do NOT put the velour cover in the freezer).

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10 Replace the pressure pad that was used during sports with the gel cushion (fig. M).

11 Attach the bandage in the same way as for sporting activities. Ensure that the gel cushion is exactly in the same position as the first pressure pad.

12 Wear the bandage for the duration prescribed by the physician.

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10 Shortly before applying the gel cushion insert it in the velour cover (fig. L).

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