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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

- REF. 460 Short tubular knee brace
- REF. 463 Short tubular knee brace with patella stabilizer
- REF. 404N Neoprene tubular knee brace with patella stabiliser
- REF. 416 Neoprene knee brace for patella stability

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. 460 - REF. 463					
	XS	S	M	L	XL	XXL
Size						
Knee center circumf. cm	27/30	31/33	34/36	37/39	40/43	44/47
Lenght cm	26	26	26	26	26	26
Colour	black REF. 460 - black and blue REF. 463					

Fits right and left

Code	REF. 404N - REF. 416					
	XS	S	M	L	XL	XXL
Size						
Knee center circumf. cm	30/33	33/36	36/39	39/43	43/48	48/53
Lenght cm	33	33	33	33	33	33
Colour	black and blue REF. 404N - black REF. 416					

Fits right and left

MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry
- Washing instructions: Remove the palm splint from its seat. Hand wash in warm water with mild soap; rinse it thoroughly. Allow to dry away from heat sources

Do not dispose of the product or any of its components into the environment.

INDICATIONS REF. 460/463/404N

- Knee pain due to bruises
- Early stage arthrosis
- REF. 460/only** First degree sprain to collateral ligaments
- Tendinitis
- REF. 463 only:**
- Chondropathy/ Chondromalacia patella
- Prevention of patella luxation during sports
- REF. 404N only**
- Chondropathy/ Chondromalacia patella
- Mild instability of the patella
- Prevention of patella luxation during sports

INDICATIONS REF. 416

- Gonarthrits
- Patellofemoral arthrosis
- Mid-lateral patella instability

CONTRAINDICATIONS

Currently no known

FEATURES AND MATERIALS

- Structure made of bi-layer orthopedic neoprene
- Lined with soft velveteen to avoid sweating and skin irritations
- Washable
- REF. 463 only:** Patella stabilizer
- REF. 404N only:** Patella pocket for U-shaped stabilizer
- REF. 416 only**
- Open kneecap with a lateral J-shaped foam stabiliser
- Velcro® fastening system
- Medial and lateral spiral stays

PUTTING ON THE APPLIANCE REF. 460/463

- Slide the tubular support onto the knee with the writing toward the top.
- Model REF. 460** (fig. A): Position the knee support so that the seams are at the side.
- Model REF. 463** (fig. B): Position the knee support so that the kneecap is perfectly centred in the opening provided.

PUTTING ON THE APPLIANCE REF. 404N

- Slide the tubular support onto the knee with the writing toward the top and position the knee support so that the kneecap is perfectly centred in the opening provided.
- If necessary, position the horseshoe stabiliser around the kneecap in the pocket provided

PUTTING ON THE APPLIANCE REF. 416

- Undo the patella straps.
- Slide the tubular support onto the knee (fig. D) and position it so that the kneecap is perfectly centred in the opening provided and the stabiliser is to the side of the kneecap (fig. E).
- Pull the patella straps and fasten them with Velcro® (fig. F-G).

N.B.: if the Velcro® is too long and pokes out, it can be cut to size.

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