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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Neoprene tubular knee brace for patella stability

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. People prone to dermatitis or allergies to products resulting from the neoprene must not use this brace. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.412				
	XS	S	M	L	XL
Size					
Circumference middle of knee cm	30/33	33/36	36/39	39/43	43/48
Length cm	33	33	33	33	33
Colour	blu				

fits right and left

MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry

 Washing instructions: Hand wash in lukewarm water with mild soap; rinse thoroughly. Leave to dry away from direct heat.

Do not dispose of the product or any of its components into the environment.

FEATURES AND MATERIALS

- It includes a standard support and a stabilising directional support with a pad that can be placed in different positions
- Multidirectional neoprene support

INDICATIONS

- Lateral-medial instability of the patella
- Femoro-patellar arthrosis
- Gonarthrits

CONTRAINDICATIONS

Currently no known.

PUTTING ON THE APPLIANCE

- 1 Position the (larger) base support with the kneecap centred in the opening provided (fig. A)
- 2 Arrange the C-shaped patella stabiliser in the desired position, and fasten to the buttons provided: there are four possible options
- 3 Pull the directional support over the part which has already been positioned so that both openings are aligned and centred over the kneecap (fig. B).

REF.412

Neoprene tubular knee brace for patella stability

object 12

**ORTHOSERVICE
RO+TEN**
Take care feel better

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**Headquarter: ORTHOSERVICE AG**
Niederlassung Deutschland: Orthoservice Deutschland GmbH
(art. 2497bis CC): **Orthoservice AG (CH)**