podifer 50

DÄMPFUNGSEINLAGE MIT QUERVERSTÄRKUNG
CUSHIONING INSOLE WITH TRANSVERSE REINFORCEMENT
SEMELLE AMORTISSEUSE À RÉNFORT TRANSVERSAL
AMORTIZIRAJUČA VKLADKA Z POPRZECZNYM WZMOCNIENIEM
AMORTIZIJOŠS ŠUPINATORS AR ŠKĒRSVIRZIENA PASTIPRINĀJUMU
PLANTARE AMMORTIZZANTE CON BARRA

REF:90050

ORTHOSERVICE
RO+TEN
 Take care feel better

Headquarter: ORTHOSERVICE AG
 Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
 Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
 info@orthoservice.com · www.orthoservice.com

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Modern ultra-thin insoles for people that spend long hours walking or standing and want more relaxed feet and enhanced comfort when walking. The special shock-absorbing material in the heel area cushions the energy produced as the foot touches the floor. An optimum combination of foam layers ensures the same effect during the rest of the stance phase. A PodoDry layer works sweat away from the foot, while the special PA and elastane covering extends the life of your insoles.
 Les caractéristiques et matériaux
 Currently no known contraindications
 1. Place the insole in the shoe and check that it is correctly positioned.
 2. Remove the original insole from the patient's shoe.
 3. Alterations may be made to the insole to suit individual requirements (it is advisable for the patient to wear the insole for a few hours before making any alterations). The foot can be shaped with scissors (use the insole that was previously in the patient's shoe for reference). The arched detailing supporting the longitudinal arch can be filed down a little at a time, as can the outer edge of the rear part of the insole.
 4. We recommend removing the insole from the shoe after wear so that it can air-dry.
 1. Place the insole in the shoe and check that it is correctly positioned.
 2. Remove the original insole from the patient's shoe.
 3. Alterations may be made to the insole to suit individual requirements (it is advisable for the patient to wear the insole for a few hours before making any alterations). The foot can be shaped with scissors (use the insole that was previously in the patient's shoe for reference). The arched detailing supporting the longitudinal arch can be filed down a little at a time, as can the outer edge of the rear part of the insole.
 4. We recommend removing the insole from the shoe after wear so that it can air-dry.

ENG. PUTTING ON THE APPLIANCE
 1. Remove the original insole from the patient's shoe.
 2. Place the insole in the shoe, check that it is correctly positioned, and ask the patient to try it on.
 3. Alterations may be made to the insole to suit individual requirements (it is advisable for the patient to wear the insole for a few hours before making any alterations). The foot can be shaped with scissors (use the insole that was previously in the patient's shoe for reference). The arched detailing supporting the longitudinal arch can be filed down a little at a time, as can the outer edge of the rear part of the insole.
 4. We recommend removing the insole from the shoe after wear so that it can air-dry.
 1. Place the insole in the shoe and check that it is correctly positioned.
 2. Remove the original insole from the patient's shoe.
 3. Alterations may be made to the insole to suit individual requirements (it is advisable for the patient to wear the insole for a few hours before making any alterations). The foot can be shaped with scissors (use the insole that was previously in the patient's shoe for reference). The arched detailing supporting the longitudinal arch can be filed down a little at a time, as can the outer edge of the rear part of the insole.
 4. We recommend removing the insole from the shoe after wear so that it can air-dry.

CONTRAINDICATIONS
 Currently no known contraindications
 CHARACTISTICS AND MATERIALS
 Top layer of fabric made with woven silver threads with natural anti-bacterial and odour-inhibiting effect.
 The PodoDry felt protects the health of the foot thanks to leading sweat away from the foot.
 The middle visco-elastic layer with form memory effect, in an optimised thickness and with damping function.
 Shock absorbing material in the heel region to cushion when setting the heel down.
 The bottom layer is made of PU and acts as a protective layer for the layers above, providing protection against wear and deformation. It also makes it easier to position the insole in footwear.
 Lower shell made of rigid PU for the stabilisation and support of the foot in pronation and supination. The presence of the transverse reinforcement supports the metatarsals and cushions the contact with the floor, in addition it strengthens the pump function of the insole, improving the venous return and, at the same time, stimulating the reflex zones for heart and lungs.
 As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

MAINTENANCE
 Wash by hand up to 30°C with neutral soap, let it dry away from heat sources. Do not dispose of the product or any of its components into the environment.
 SAFETY PRECAUTIONS
 We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or warts. It is advisable to not over-tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.
 WARNINGS
 It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibility of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, warts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

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Modern ultra-thin insoles for people that spend long hours walking or standing and want more relaxed feet and enhanced comfort when walking. The special shock-absorbing material in the heel area cushions the energy produced as the foot touches the floor. An optimum combination of foam layers ensures the same effect during the rest of the stance phase. A PodDry layer wicks sweat away from the foot, while the special PA and elastane covering extends the life of your insoles.

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