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SINGLE PATIENT  
MULTIPLE USE

 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

## Low/high elastic orthosis with a rigid thermoformable pad to stabilize the lumbar rachis and relieve pressure

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademarc of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin.

If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

### WARNINGS


It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Item	REF. 1770H	REF. 1770S	REF. 1870H	REF. 1870S		
Size	<b>XS</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
Circ. waist cm	50/60	60/75	75/90	90/105	105/120	120/135
Front height cm	20	20	20	20	20	20
Back height cm	26	26	26	26	26	26
Colour	REF. 1870H - REF. 1770H = light blue / REF. 1870S - REF. 1770S = black					

Item	REF. 1772H	REF. 1772S	REF. 1872H	REF. 1872S		
Size	<b>XS</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
Circ. waist cm	50/60	60/75	75/90	90/105	105/120	120/135
Front height cm	22	22	22	22	22	22
Back height cm	32	32	32	32	32	32
Colour	REF. 1772H - REF. 1872H = light blue / REF. 1772S - REF. 1872S = black					

### MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry
-  Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources. Dispose of this product and its parts responsibly

### INDICATIONS

#### Linear 70/72

- Lumbalgia
- Paravertebral muscular contractures
- Linear70 only: Mild traumas of the lumbosacral vertebral column
- Linear72 only: Mild traumas of the lumbar vertebral column

#### Linearplus 70/72

- Lumbago
- Lumbar sciatica and lumbar cruralgia
- Spondyloarthritis and lumbar discopathy
- Slight traumas to the lumbosacral column and paravertebral lumbar contractures
- Postoperative recovery
- Linearplus72 only: Results of vertebral compressions of the lumbar spine
- Linearplus72 only: Treatment of transverse process fractures

### KONTRAINDICATIONS

Currently no known

### CHARACTERISTICS AND MATERIALS

#### Linear 70/72

- Structure made from ribbed cotton and elastomer fabric
- Front Velcro® fastening system
- Linear reinforcement on the back
- Adapts to the pathology of the patient

#### Linearplus 70/72

- Structure made from ribbed cotton and elastomer fabric
- Front Velcro® fastening system
- Crossed tensioning straps at the back to increase the compression on the lumbar region
- Adapts to the pathology of the patient

### PUTTING ON THE APPLIANCE

#### Linear 70/72

- Open the orthosis and position with the label facing upwards.
- Hold the orthosis for fi tting and fi x well using the Velcro® closure. Ensure that the closure is located exactly in the middle (fi g. A).
- In case of prescription of a lumbar pelotte, fix it on the back at the Velcro® closure prior to fitting the orthosis.

#### Linearplus 70/72

- Open the orthosis and temporarily fi x x the Velcro® straps at the lateral straps (fi g. B)
- Open the orthosis and position with the label facing upwards.
- Hold the orthosis for fi tting and fi x well using the Velcro® closure. Ensure that the closure is located exactly in the middle (fi g. C).
- Tighten the orthosis by fixing the Velcro® straps on the front Velcro® strips. ONLY high model: First pull tight the lower and then the upper Velcro® straps (fi g. D).
- In case of prescription of a lumbar pelotte, fix it on the back at the Velcro® closure prior to fitting the orthosis.

### NOTE:

To ensure optimal stabilisation, mould the rigid stays provided and slide them into the pockets provided in the back.

#### Attaching the rigid pad:

- if necessary, heat and mould the plastic part to fit the shape of the patient's body
- place the padding on the plastic part (concave side-without any writing-fi g.E)
- apply the adhesive Velcro® fasteners centrally on the opposite side to the padding (fi g.F)
- attach the pad to the central Velcro® to the inside of the back of the orthosis

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.