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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Arm and shoulder immobilizer

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS


The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Item	REF.8410	
Size	Regular	Large
Waist circumference in cm	< 120	> 120
Colour	black	

Fits right and left

MAINTENANCE

 Do not bleach No chemical cleaning
 Do not iron Do not tumble-dry
 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Stable humerus fractures
- Gleno-humeral luxation
- Distorsions of shoulder
- Post surgery humerus fractures and gleno-humeral luxation treatment
- Post surgery after tenorrhaphy of long head of humeral biceps

CONTRAINDICATIONS

Currently no known

FEATURES AND MATERIALS

- Soft structure with external gauze nylon to match Velcro® fastening system
- Shoulder strap with snap closure

APPLICATION

- Pass the waist belt (A) around the waist and use the Velcro® closure to tension firmly and close.
- Bring the arm cuff (B) that is sewn to the waist belt up to the height of the affected arm and then secure by closing the appropriate Velcro® closure.
- After inserting the arm, secure the lower arm support brace (C) to the waist belt (A).
- On the rear, secure the closure strap that encircles the shoulder (D) to the waist belt (A) using the Velcro® closure. Pass the strap (D) over the shoulder opposite of the affected one and secure to the lower arm support brace with the clip buckle closure (E).
- Tighten the strap (D) appropriately. In doing so, ensure that the protective pad for the neck (F) is positioned correctly.
- Fix the upper part of the lower arm support brace (C) to the waist belt (A) using the Velcro® securing strap (G).

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